

# Maximum Procedure Flexibility

## Secour® Acrylic Resin



Parallax Secour Acrylic Resin with Monomer

**With Secour Acrylic Resin, save time, procedure costs and improve patient satisfaction.**

### Unique Features:

- Easy to use with bone cement opacifiers including TRACERS® – just add opacifier, monomer and mix together
- A 12-15 minute working time increases procedure flexibility
- Pre-measured single dose package of 12g of polymer and 7.2ml ampoule of Monomer ensures accurate and consistent mixing every time
- Indicated for use with TRACERS, a proprietary barium sulfate opacifier. TRACERS provides optimum visualization for the placement of bone cement in a vertebral body and does not alter cement working characteristics or compression strength

*Simplifying Procedures. Satisfying Patients.*

# Secour® Acrylic Resin



Pre-measured single dose packaging

*Secour is part of the line of leading acrylic bone cement products and opacifiers from ArthroCare Spine.*

#### These products include:

- TRACERS®: barium sulfate bone cement opacifier
- Parallax Acrylic Resin with TRACERS: Secour plus TRACERS opacifiers in one easy to use and mix package. FDA indicated for vertebroplasty and kyphoplasty
- Parallax Acrylic Resin with TRACERS-Ta: Secour plus our exclusive formula of Tantalum opacifiers and barium sulfate in one easy to use and mix package

**Premeasured polymer in a single dose package ensures accurate, consistent mixing, repeatable performance and less material waste**

- Packaged as 12g of powder in a mixing vial and 7.2ml of liquid monomer in an ampule
- Secour is a resin PMMA (poly-methylmethacrylate) cement
- Meets ASTM F451 and ISO 5833 requirements for compressive strength, implants, intrusion depth and maximum exothermic temperature for acrylic resin implants

#### ORDERING INFORMATION

PRODUCT DESCRIPTION	Secour Acrylic Resin PMMA
CATALOG NUMBER	SAR-901
CONTENTS	12g of polymer powder in a mixing vial 7.2ml of liquid monomer in an ampule

*See package insert for indications, warnings, precautions and instructions for use.*

To order or for more information, please contact your ArthroCare Spine representative today.



**ArthroCare SPINE**  
680 Vaqueros Avenue  
Sunnyvale, CA 94085-3523  
toll-free **800-797-6520**  
phone 408-736-0224  
fax 408-736-0226  
[www.parallax-medical.com](http://www.parallax-medical.com)

**ArthroCare Europe AB**  
Baggensgatan 25  
111 31 Stockholm  
Sweden  
phone +46 8 546 172 00  
fax +46 8 546 172 39  
[info@arthrocare.se](mailto:info@arthrocare.se)

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

© 2005 ArthroCare Corporation. ArthroCare, Parallax, Secour and TRACERS are trademarks and/or registered trademarks of ArthroCare Corporation. US Patents 6019776, 6033411, 6231615, 6309420, 6348055, 6383190 and patents pending.